



Instituto Tecnológico de Estudios Superiores de Monterrey
Campus Monterrey

Procesamiento de imágenes médicas para el diagnóstico (Gpo 200)

MRI Activity

Profesor:

Dr. José Gerardo Tamez Peña

Equipo 3:

Ian Salvador Becerra Sánchez A01274634

Elías Marcelo Martínez Muñozcano A00829233

Karime Dagdug Chacón A00827240

Eugenio Rubio Cantú A00829391

Manuel Verduzco Valenzuela A00827326

Monserrat Adely Izaguirre A01412255

MRI activity.

Open the simulator in <https://phet.colorado.edu/en/simulations/mri>

1. NMR

Place yourself in the tab *Simplified NMR*. Remember that the Larmour relationship relates the magnetic field to the resonant frequency:

$$\nu_L = \frac{\gamma}{2\pi} B_0$$

where ν_L is the resonant frequency, $\frac{\gamma}{2\pi}$ is the gyromagnetic ratio and B_0 is the magnetic field.

Table1. Gyromagnetic ratios for different nuclei.

Nuclei	Gyromagnetic ratio	Nuclei	Gyromagnetic ratio
¹ H	42,58	⁶⁵ Cu	12,09
⁷ Li	16,55	⁷⁵ As	7,291
⁹ Be	5,984	⁷⁷ Se	8,118
¹¹ B	13,66	⁸¹ Br	11,50
¹³ C	10,71	⁸⁷ Rb	13,93
¹⁵ N	4,314	⁹³ Nb	10,41
¹⁷ O	5,772	¹¹⁷ Sn	15,17
¹⁹ F	40.05	¹²¹ Sb	10,19
²³ Na	11,42	¹²⁷ I	8,518
²⁷ Al	11,09	¹³³ Cs	5,584
²⁹ Si	8,458	¹⁹⁵ Pt	9,153
³¹ P	17,24	¹⁹⁹ Hg	7,590
³⁵ Cl	4,172	²⁰³ Tl	24,33
⁵¹ V	11,19	²⁰⁷ Pb	8,907
⁵⁵ Mn	10,50	²⁰⁹ Bi	6,841
⁵⁹ Co	10,05		

Use the Larmour relationship and the gyromagnetic ratios of various nuclei shown in table 1 to complete table 2. Check your results in the simulation by setting the appropriate frequencies and main magnet field, take a screenshot of the nuclei emitting energy to include in the report. Try to find the last nuclei (????) by playing with the simulation and register the frequency at two different magnetic fields.

Table 2. Different settings to achieve energy emission.

Nuclei	Magnetic Field	Resonant Frequency	Magnetic Field	Resonant Frequency
Hydrogen	0.75	31.935	2.5	106.45
Nitrogen	2.5	10.785	-	-
Sodium	1.5	17.3	2.75	31.405
Carbon-13	1.75	18.7425	2.5	26.775
Oxygen	2	11.54	3.0	17.316
Sulfur	3.0	10	-	-
????				

2. MRI

Move to the *Simplified MRI* tab

- Set the *main magnet field* to 1.0 Tesla, leave the *gradient magnets* in zeros, activate only *show head*, and *show magnetic field* (be sure that *show atomic nuclei* is deactivated), set the frequency in 43 MHz. Finally set the *power* to 50% and observe the flow and distribution of the emissions. After a while observing the emissions, click on *add tumor*, wait for around 7 seconds so the flow distribution stabilizes, look at how the emission changed and try to guess where the tumor is located.

Explain how the emission allowed you to find the correct location: The tumor is located on the top right side of the head. We can determine this because of how the emission changes when the tumor is added. We can see that there is a higher wave concentration, this is what helps us find the location of the tumor.

- Play with the main magnet field, frequency, and gradient magnets (both, horizontal and vertical) to try to obtain an emission focused mainly in the zone of the tumor (register your best guess, it doesn't need to be perfect). Answer the following questions.

Best guess: main magnet: 0.5 horizontal gradient: 0.08 vertical gradient: 0.01 frequency: 20
 What happens when the horizontal gradient increases its magnitude? How does it affect the emissions? There is a greater concentration of emission in the x axis, depending of the value of the gradient.

What about vertical gradient? It affects the concentration of emissions on the y axis.